

Nava G, Adams JA, Bassuk J, Wu D, Kurlansky P, Lamas G. Echocardiographic comparison of cardiopulmonary resuscitation (CPR) using periodic acceleration (pGz) versus chest compression. Resuscitation 2005;66:91-7

Divisions of Cardiology, Neonatology, Department of Research, Mount Sinai Medical Center, Miami Heart Research Institute, 4300 Alton Road, Miami Beach, FL 33140, USA.

OBJECTIVE: This investigation compared the effects of conventional cardiopulmonary resuscitation (CPR) using an automated Thumper chest compression device to periodic acceleration CPR (pGz-CPR) on early post-resuscitation ventricular function assessed by echocardiography, in an adult pig model of CPR. **BACKGROUND:** Whole body periodic acceleration along the spinal axis (pGz) is a new method of cardiopulmonary resuscitation (CPR). Biomechanical forces and biochemical release produced by pGz impart ventilation and increase blood flow. Our laboratory has reported normal neurological and cardiovascular function 48 h after return of spontaneous circulation in animals that have undergone 22 min of pGz-CPR. **METHODS:** Ventricular fibrillation (VF) was induced in 16 animals (25-35 kg). After 3 min of non-interventional period, the animals were randomized to receive either pGz-CPR or Thumper-CPR for 15 min. After 18 min of VF, a single dose of vasopressin and bicarbonate were administered and defibrillation attempted. An echocardiogram was performed at baseline and serially for 6h. Ejection fraction (EF), fractional shortening (FS) and wall motion were assessed by 2D and M-mode echocardiography. **RESULTS:** Return of spontaneous circulation to 360 min occurred in 5/8 (62%) of the animals receiving Thumper-CPR and in 7/8 (88%) receiving pGz-CPR. FS and EF were impaired after CPR, but pGz-CPR animals had less impairment than Thumper-CPR animals. Further, wall motion score index (WMSI) was more impaired after Thumper-CPR and remained as such even 6h post-CPR. **CONCLUSION:** pGz holds promise as a new method for CPR with better left ventricular (LV) function post-CPR than the more traditional chest compression method.